MONARK
894 E

OBS!

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NOTE!

Use of the product may involve considerable physical stress. It is therefore recommended that people who are not accustomed to cardiovascular exercise or who do not feel completely healthy, should consult a physician for advice.

When using a chest belt, keep in mind that electromagnetic waves can interfere with the pulse signal. Mobile phones and the like should therefore not be used near the bike when using a chest belt.


For service, maintenance and spare parts list, refer to our service manual. You can also contact your dealer.

Max user weight 125 kg.





Before using your new bike, please read the user information carefully and save it for future reference.

NOTE!


If you want to perform aerobic or anaerobic tests using this product, you should have relevant knowledge. We therefore refer to the relevant literature that describes the test procedure and analysis of measured results.

If the product or its packaging is marked with this symbol, it should not be handled as household waste. Instead, it shall be provided to a suitable collection point for the recycling of electrical and electronic equipment.


5)


NOTE!

The stem must be inserted at least 8 cm in the frame. This is marked with "MAX" on the stem.

6)

2 x 1.5V, R6(AA)

7)


1x R

1x L


1x

NOTE!

When the bike is new, or if you replace the pedals, they must be tightened after 5 hours of use.

8)

DISPLAY

Press any button or start pedalling to start up the display. Then the display start counting the values for rpm, time, speed and distance and pulse (HR) (if chest belt is used).

In addition, the braking force (kp = weight basket + weights in kg) can be set which gives a correct indication of calorie consumption (CAL) and power (WATT). The effect is dependent on the rpm, which means that the wattage can be fine-tuned by increasing alt. reduce rpm.

At HR, a heart symbol lights up, which means that the display searches for a pulse signal. This is disconnected after 30 seconds unless a chest belt is found. You can start the search again by pressing the RESET button. When a pulse signal is found, the heart rate is shown in the display.

Time begins to count automatically as the pedals are pedalled. The values for Time (TIME), distance (DISTANCE) and calories (CAL) can be reset by pressing the RESET button for more than two seconds.

In order for the display of calories and watts to show correctly, the kp (work load) must be set to the load value corresponding to the weight basket (weights + weight of the basket 1 kg). The rubber plates are included in the calibrated weight for 1 kg weight basket.

Example:

The workload is 3 kg (weight basket 1 kg + 2 x 1 kg weight). Press the kp button. The lower display window now shows the kp value. Increase in increments of 0.1 with the kp button, alt. decrease with the RESET button until the value matches the total weight of the weight basket, in this case 3.0. Then press the CAL / WATT button and the cycle's braking power in watts is displayed. At further press on the CAL / WATT-button the display toggles between calories and watts. The wattage displayed will depend on the rpm and the wattage can therefore be fine-tuned by increasing or decreasing the rpm.

Calorie consumption is recalculated continuously.

NOTE!

Do not expose the electronic display to direct sunlight or extremely high temperatures. Do not use solvents when cleaning, only a dry cloth.



Fig:
Electronic display

Change units

On the back of the display is a switch for displaying units in km / h and km or display in mph and miles. Then, remove the batteries (without using tools) and reinsert them.

1 = km/h and km

ON = mph and miles

WORKLOAD ADJUSTMENT

Changing the work load is done either by another pedalling speed, or by increasing or decreasing the tension of the brake belt / cord against the brake wheel by putting weights in the weight basket. Weights are available in 1 kg, 0.5 kg and 0.1 kg. This makes it possible to vary the workload from 1 kp up to maximum 11 kp in steps of 0.1 kp. Braking power is expressed in kp where mass 1 kg gives the braking force 1 kp.

If you do not want any work load at all, the weight basket can be hung in its upper stop position. The weight basket is released - work load is applied - by pressing the push button located on the end of the handlebar. See *Fig: Workload adjustment*.


Fig:
Workload adjustment
1) Release button, release basket
2) LED

NOTE!

1 kg is the lowest work load that can be set as this is the weight of the basket itself. A weight basket that only weighs 0.5 kg is available as an accessory.

HEART RATE

The chest belt belt is supplied as standard.

In order for the heart rate measurement to work, it is necessary that the pulse belt is correctly positioned, see *fig: Placement chest belt*. Moisten the electrodes before use, see *Fig: Moistening the electrodes*.

To make contact with the HR receiver on the bike, the distance should not be more than 100 cm. At the first usage it is important to first pair the chest belt with the display by standing close to it so the display can receive the pulse signal (max. 60 cm).


Fig:
Placement chest belt


Fig:
1) Moistening the electrodes

NOTE!


Electromagnetic waves may disturb heart rate function. Cellular phones are not allowed to be used near the bike during test.

CONNECTION TO PC

To connect a PC to the bike, use a 0-modem cable (RS232) with 9-pin D-sub connector (female) at both ends. If no RS232 Serial port is available on the computer use a USB serial converter to connect to an USB port on the computer instead.

Before installing the Monark Software for your ergometer, you must take the following steps.

1. Locate the USB adaptor.
2. Inside the USB adaptor packaging, there is a mini-CD.
3. Insert the mini-CD into the CD-ROM drive and install the driver software. If there is no CD drive on your computer, driver software is available for download from the website (<http://www.vscom.de/USB-CD>).
4. Finalize the driver software installation by inserting the USB adaptor.
5. Install the Monark Software disk.
6. Connect the USB Adaptor to the Serial Cable and proceed with testing.

ANAEROBIC TESTING

The ergometer contains a computer that allows it to be connected to an external PC via a cable that is delivered with the ergometer. The software - Monark ATS Software - you can download from our website, www.monarkexercise.se. It makes it possible to easily run a variety of different types of anaerobic tests, eg. Wingate tests. The duration of the anaerobic tests can easily be set from five seconds up to 300 seconds. For a more detailed description of the software, see separate program manual "Monark Anaerobic Test Software - User Manual".

To carry out anaerobic tests must cycle computer equipped with power from a wall outlet. Connect the supplied power adaptor into an appropriate wall outlet (230V). The circular connector on the cable from the AC adapter is connected to the corresponding connector on the ergometer (under the instrument cover on the right side). The PC cable connects to the data output on the underside. See *Fig: Brake device*. The LED on the top of the instrument panel indicates both that power is on the unit and that the weight basket is in its upper stop position / neutral position. If the weight basket produces resistance it must be moved up in its upper locked position where it does not give any workload. Then the LED is lit - with a delay of about five seconds - indicating that the computer and cycle is in a ready position for a test.

Fig:
Brake device
1) Release button
2) Stop
3) Cord bracket
4) Tension center
5) Return strap
6) Weight basket
7) DC connector
8) PC serial port
9) Pin on tension center
10) Tip of the magnet arm

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